



CL007N English

CL007N Roller Mixer

Our new generation Roller Mixer with a steeper angle and larger rolls enables an even better mixing at continuous speed.

Choose **CYAN CL007N** and get the following advantages:

- Suitable for all conventional tubes and cylindrical bottles
- Safe and quiet instrument
- Ellipsoidal orbit
- Smooth rotation avoids foaming of the samples
- Ideal for maintaining blood sample quality prior to analysis
- Rocking and rolling movement

Technical data:

- Speed: 40 rpm
- Roll dimension: 5 x 32 cm (Ø x length)
- 5 rolls
- Overall dimension: 45,5 x 28,3 x 9,7 cm (W x D x H)
- Gross weight: 7,4 kg
- Power supply: 110 V - 60 Hz (CL007N-110)
or 220 V - 50/60 Hz (CL007N-220)

Order Code:

- CL007N-110 [Roller Mixer](#) (110 V - 60 Hz)
- CL007N-220 [Roller Mixer](#) (220 V - 50/60 Hz)

Roller Mixer



Contact: e-mail: cypress@diagnostics.be
tel: + 32 (0)15 67 67 68

Headquarters: Cypress Diagnostics
Langdorpsesteenweg 160
3201 Langdorp
Belgium

Factory: Nijverheidsstraat 8
2235 Hulshout
Belgium

URL: www.diagnostics.be

Your Distributor

Roller Mixer EN 2018-01