



CL016 English

CL016 Mini Vortex

Suitable for mixing the contents of specimen tubes in the laboratory.
The instrument operates automatically upon pressing on the head of the vortex (press start function).

Choose **CYAN CL016 and get the following advantages:**

- Automatic press start function
- Variable speed control
- Compact, robust, and durable
- Vibration free
- Fast
- Quiet
- High-quality motor for perfect performance

Technical data:

- Silicon cup head
- Orbit of 4,5 mm
- Max. speed: 3200 rpm
- Max. volume: 50 ml
- Max. weight: 60 g
- Head cup dimensions: 3,9 cm (Ø)
- Overall dimensions: 5,5 x 11,3 x 8,1 cm (W x D x H)
- Gross weight: 0,90 kg
- Power supply: 110/220 V - 50/60 Hz

Order code:

- CL016: **Mini Vortex** (110/220 V - 50/60 Hz)

Mini Vortex



Contact: e-mail: cypress@diagnostics.be
tel: + 32 (0)15 67 67 68

Headquarters: Cypress Diagnostics
Langdorpsesteenweg 160
3201 Langdorp
Belgium

Factory: Nijverheidsstraat 8
2235 Hulshout
Belgium

URL: www.diagnostics.be

Your Distributor

Mini Vortex EN 2018-01