



CL010 English

### CL010 Water bath

Choose **CYAN CL010** and get the following advantages:

- Low electricity consumption
- Very stable temperature
- Bath with rounded corners for easier cleaning

**Technical Data:**

- Unstirred (thermal diffusion)
- Working range: ambient temperature up to 100 °C
- Capacity: 10 liters
- PID temperature control
- Digital temperature setting
- Digital temperature readout
- Stability: +/- 0,2 °C
- Inner tank dimensions: 24 x 30 x 15 cm (W x D x H)
- Overall dimensions: 29 x 40 x 24 cm (W x D x H)
- Gross weight: 8,5 kg
- Overheating protection
- Power: 650 W
- Power supply: 110 V – 60 Hz (CL010-110) or 220 V – 50/60 Hz (CL010-220)

**Accessories included:**

- Flat cover for water bath

**Order code:**

- CL010-110: [Water bath](#) (110 V – 60 Hz)
- CL010-220: [Water bath](#) (220 V – 50/60 Hz)

**Auxiliary products:**

- CL010-R Test tube rack 50 places



# Water Bath



Contact: e-mail: [cypress@diagnostics.be](mailto:cypress@diagnostics.be)  
tel: + 32 (0)15 67 67 68

Headquarters: Cypress Diagnostics  
Langdorpsesteenweg 160  
3201 Langdorp  
Belgium

Factory: Nijverheidsstraat 8  
2235 Hulshout  
Belgium

URL: [www.diagnostics.be](http://www.diagnostics.be)

Your Distributor

Waterbath EN 2018-01